

HearingLoop NW

P. Spencer Norby
HearingLoop NW
2522 N. Proctor #519
Tacoma, WA 98406
253.686.6684
psnorby@hearingloopnw.com



www.hearingloopnw.com

Loop Instructions:

A hearing loop has been installed in our Meeting Room


Hearing aid, cochlear implant and BAHA wearers can switch their hearing devices to “T-Coil” mode and hear the sound directly in their hearing devices.

If you don’t know what that is, most hearing appliances have a coil inside them that picks up the hearing loop directly. The hearing loop gives off a vibrating magnetic signal. The t-coil picks up this magnetic signal and turns it back into sound.

If you don’t know if you have a “T-Coil”, cycle through your program button on your hearing device and see if the sound changes so it seems to be right in front of you. The program button is usually small and does several different things. It will probably beep each time it is pressed. Some hearing devices have a remote control for them that has a “T-Coil” in the remote. The next page has pictures of program buttons.

If your hearing device does not seem to sound any different when you cycle through, you will have to check with your audiologist. You may have a “T-Coil” that just isn’t activated.

You may have to ask about having your device fitted with a “T-Coil”. If that can’t be done, explain to your audiologist that you need a “T-Coil” and you want them to figure out how to get you one.



VHP-220T

